

	MacPractice 7.3			MacPractice Gen 10		
	El Capitan (10.11)	Sierra (10.12)	High Sierra (10.13)	Sierra (10.12)	High Sierra (10.13)	Mojave (10.14)
3rd Party Applications						
Microsoft Office 2011	Yes	Yes	No	Yes	No	No
Microsoft Office 2016	Yes	Yes	Yes	Yes	Yes	Yes
QuickBooks	Yes	Yes (QuickBooks 2016)	Yes (QuickBooks 2016)	Yes (QuickBooks 2016)	Yes (QuickBooks 2016)	Yes (QuickBooks 2016)
Quicken	Yes	Yes (Quicken 2016+)	Yes (Quicken 2016+)	Yes (Quicken 2016+)	Yes (Quicken 2016+)	Yes (Quicken 2016+)
Parallels	Yes (Parallels 10-12, Windows XP - 10)	Yes (Parallels 11-12, Windows XP - 10)	Yes (Parallels 13, Windows 8 - 10)	Yes (Parallels 11-12, Windows XP - 10)	Yes (Parallels 13, Windows 8 - 10)	Yes (Parallels 13, Windows 8 - 10)
VMWare Fusion	Yes (Fusion 8.0-8.5, Windows XP - 10)	Yes (Fusion 8.5, Windows XP - 10)	Yes (Fusion 8.5, Windows XP - 10)	Yes (Fusion 8.5, Windows XP - 10)	Yes (Fusion 8.5, Windows XP - 10)	Yes (Fusion 11+, Windows XP - 10)
MacPractice Services						
AutoRemind	Yes	Yes	Yes	Yes	Yes	Yes
CDA ITrans	Yes	Yes	Yes	Yes	Yes	Yes
Change Healthcare	Yes	Yes	Yes	Yes	Yes	Yes
eStatements (Change Healthcare)	Yes	Yes	Yes	Yes	Yes	Yes
Integrated Network Faxing	Yes	Yes	Yes	Yes	Yes	Yes
Online Patient Registration	Yes	Yes	Yes	Yes	Yes	Yes
Online Patient Scheduling	Yes	Yes	Yes	Yes	Yes	Yes
Patient Portal	Yes	Yes	Yes	Yes	Yes	Yes
Reputation Marketing	Yes	Yes	Yes	Yes	Yes	Yes
Secure Messaging	Yes	Yes	Yes	Yes	Yes	Yes
Transworld Collections	Yes	Yes	Yes	Yes	Yes	Yes
AuthPay X	Yes	Yes	Yes	Yes	Yes	Yes
TSYS (TransFirst Credit Card Interface)	Yes	Yes	Yes	Yes	Yes	Yes
DollyDrive Internet Backup	Yes	Yes	Yes	Yes	Yes	Yes
X-Ray Native Interface						
These native interfaces send images directly to the Digital Radiography ability in MacPractice						
Acteon SOPIX2 USB Sensor (1, 2)	Yes (Only on Macs from 2012 and newer)	Yes (Only on Macs from 2012 and newer)	Yes	Yes (Only on Macs from 2012 and newer)	Yes	Yes
Acteon PSPIX Phosphor Plate (0, 1, 2, 3, 4)	Yes (Only on Macs from 2012 and newer)	Yes (Only on Macs from 2012 and newer)	Yes	Yes (Only on Macs from 2012 and newer)	Yes	Yes
DEXIS Mac*	Yes	Contact Vendor	Contact Vendor	Contact Vendor	Contact Vendor	Contact Vendor
MyRay X-pod with Sensor (1, 2)	Yes (Bluetooth)	Yes (Bluetooth)	Yes (Bluetooth)	Yes (Bluetooth)	Yes (Bluetooth)	Yes (Bluetooth)
Planmecca ProSensor (0, 1, 2)	Yes	Yes	Contact Vendor	Yes	Yes	Contact Vendor
Planmecca Romexis	Yes	Yes	Contact Vendor	Yes	Yes	Contact Vendor
QuickRay USB Sensor (1, 2)	Yes	Yes	Yes	Yes	Yes	Yes
SuniRay/SuniRay 2 Sensor (1, 2)	Yes	Yes	Yes	Yes	Yes	Yes
* DEXIS Mac stores images in a separate local server and not in MacPractice DR						
X-Ray Native Interface (Pan/Ceph) Cone Beam						
MyRay Hyperion	Yes (iPad to OS X)	Yes (iPad to OS X)	Yes (iPad to OS X)	Yes (iPad to OS X)	Yes (iPad to OS X)	Yes (iPad to OS X)
Planmecca ProOne and ProMax	Yes	Yes	Contact Vendor	Yes	Contact Vendor	Not Yet
X-Ray Bridge (Windows)						
Please verify your computer or VM Software can run the appropriate version of Windows that your Digital Radiography Software Requires.						
Apteryx/CliniView (Most Sensors, Pan, Ceph)	Yes	Yes	Yes	Yes	Yes	Yes
CADI	Yes	Yes	Yes	Yes	Yes	Yes
DEXIS	Yes	Yes	Yes	Yes	Yes	Yes
Gendex Sensor / Pan / Ceph (VixWin)	Yes	Yes	Yes	Yes	Yes	Yes
Instrumentarium Pan / Ceph	Yes	Yes	Yes	Yes	Yes	Yes
Lightyear	Yes	Yes	Yes	Yes	Yes	Yes
OpTime (Phosphor Plate)	Yes	Yes	Yes	Yes	Yes	Yes
ProfSuni	Yes	Yes	Yes	Yes	Yes	Yes
ScanX / Visix (Phosphor Plate)	Yes	Yes	Yes	Yes	Yes	Yes
Schick (Pre-Sirona)	Yes	Yes	Yes	Yes	Yes	Yes
Sirona Sidexis 4	Yes	Yes	Yes	Yes	Yes	Yes
VixWin / Gendex	Yes	Yes	Yes	Yes	Yes	Yes
Imaging Software Bridge (Mac)						
Planmecca Romexis	Yes	Yes	Yes	Yes	Yes	Contact Vendor
RadioVision	Yes	Yes	Yes	Yes	Yes	Yes
Perio Charting Data Entry - Foot Pedal						
Dental R.A.T.	Yes	Yes	Yes	Yes	Yes	Yes
Camera (Plug and Play)						
DentSmart Dental Camera	Yes	Yes	Yes	Yes	Yes	Yes
iCam	Yes	Yes	Yes	Yes	Yes	Testing
ImageMaster Billet USB Intraoral Camera	Yes	Yes	Yes	Yes	Yes	Yes
ImageMaster USB	Yes	Yes	Yes	Yes	Yes	Yes
MouthWatch	Yes	Yes	Yes	Yes	Yes	Yes
Camera (With Interface Software)						
These cameras can be used with MacPractice by using interface software (listed below) and saving images to an import folder						
Digital Doc Iris 1000 (Standard)	No	No	No	No	No	No
Digital Doc Iris Mac Version	Yes	Yes	Yes	Yes	Yes	Some compatibility issues
Discovery Cam	Yes	Yes	Yes	Yes	Yes	Some compatibility issues
Discovery Cam 2	Yes	Yes	Yes	Yes	Yes	Some compatibility issues

DocPortMacro USB	Yes	Yes	Yes	Yes	Yes	Some compatibility issues
DR Cam	Yes	Yes	Yes	Yes	Yes	Some compatibility issues
SoproCare	Yes	Yes	Yes	Yes	Yes	Yes
SoproLife	Yes	Yes	Yes	Yes	Yes	Yes
Sopro 617	Yes	Yes	Yes	Yes	Yes	Yes
Sopro 717	Yes	Yes	Yes	Yes	Yes	Yes
SUNI SuniCam II USB	Yes	Yes	Yes	Yes	Yes	Yes
Camera Interface Software						
ProView USB Software	Yes	Yes	Some compatibility issues	Yes	Some compatibility issues	Some compatibility issues
SoproCapture	Yes	Yes	Yes	Yes	Yes	Yes
VideoGlide	Yes	Yes	Some compatibility issues	Yes	Some compatibility issues	Some compatibility issues